

3PLY SURGICAL MASK

PN: 020001 / MANUFACTURER PN: TNK-WKA

EN 14683:2019 Type IIR



Tested and validated by



SPEC

Dimension	17.5cm ± 0.5cm x 9.5cm ± 0.5cm
Outer layer:	nonwoven PP 25gsm ± 10%, blue
Filter layer:	melt-blown 25gsm ± 10%, white, BFE 98%
Inner layer:	nonwoven PP 25gsm ± 10%, white

Packaging

Plastic bag	10PCS, 23 x 13cm
Box	50PCS, 5 packs, 22.5 x 13.2 x 8.8cm
Carton	3000PCS, 60 boxes, 63 x 50 x 49cm

