

FFP2 PROTECTIVE MASK

PN: 010103 / MANUFACTURER PN: HZ-KN95

Economic foldable FFP2




CE Certificate	
Module B	2163-PPE-623
Module C2	2163-PPE-623/01

SPEC

Standard	EN 149:2001+A1:2009
Type	FFP2 Non-reusable
Filter efficiency:	≥ 94%

Material

1st layer	non-woven PP 50gsm, white
2nd layer	melt-blown 25gsm, white, FE 95%
3rd layer:	Antibacterial static cotton, 40gsm
4th layer:	non-woven PP 30gsm, white



Packaging

Plastic bag	5 PCS
Box	10 packs
Carton	600PCS, 12boxes, 66 x 43 x 41cm

