

certificate
no. IFA 1601104
dated 2016-10-21

Translation In any case, the German original shall prevail.

European notified body
Identification number: 0121

EC Type-Examination Certificate

Name and address of the holder of the certificate (customer):
Wuxi Sunli Industrial (PSP) Co., Ltd.
BLK.B(11B), De Yuan Building
78 Wu Ai Road, Wuxi, Jiangsu
PEOPLE'S REPUBLIC OF CHINA

Product designation: **TE YIN JTY B-11F FFP3 NR D EN149:2001+A1:2009**

Type: Particle filtering half mask with exhalation valve

Testing based on: DIN EN 149:2009

Test report: 201622050/2120 of 2016-10-21 - IFA, Sankt Augustin

Further details: For protection against non highly volatile liquid and solid particles.

Class of device: FFP3 NR D

Maximum application time or maximum particle exposure one single shift.


This document is valid in connection with the technical file and the instructions for use in German language, both provided with the IFA endorsement 1601104 of 2016-10-21.

Manufacturer:
Suzhou Teyin Non-woven Fabrics Co., Ltd.
No. 2163, Qiongling Road, Zangshu Street
Mudu Town, Wuzhong District Suzhou, Jiangsu
People's Republic of China

The product tested complies with the essential requirements of Annex II of the directive 89/686/EEC (**Personal Protective Equipment**), amended by the directives 93/68/EEC, 93/95/EEC and 96/58/EC.

The present certificate is valid until: **2021-10-20**

Further provisions concerning the validity, the extension of the validity and other conditions are laid down in the Rules of Procedure for Testing and Certification.



Dr. rer. nat. Peter Paszkiewicz
Head of testing and certification body



Dipl.-Chem. Christoph Jochen
Certification officer