

# FFP3 PROTECTIVE MASK WITH VALVE

PN: 016615 / MANUFACTURER PN: JTY B-11F

CE certified by DEKRA and IFA



### CE Certificate

Module B	IFA1601104
Module C2	ZQ/A256/19

Tested and validated by




### SPEC

Standard	EN 149:2001+A1:2009
Type	FFP3 Non-reusable
Filter efficiency:	≥ 99%

### Material

Outer layer	Shaped Cotton
Mid layer	Soft Cotton
Mid layer:	P3 level Melt Blown Filter, 90GSM
Inner layer:	PP non-woven



### Packaging

Box	10 PCS
Carton	20 boxes

