


# CPE GOWN

**PN: 60002 / MANUFACTURER PN: CN615**

## Blue-Thumb-loop cuff, Waist ties

The CPE Gown is ideal for hospital and clean area usage, light and permeable, showing low resistance to water vapour permeability, especially designed to allow comfort during long duration surgery.

Apron style, long sleeves, back empty, thumbs up/elasticated wrists, 2 ties.



Size	Universal size
A	104.1 ± 3.8
B	111.8 ± 2.5
C	83.8 ± 2.5
D	208.3 ± 5.1
E	31.8 ± 1.3
F	88.9 ± 1.3
G	116.8 ± 2.5
H	68.6 ± 2.5
I	30.5 ± 1.3
J	7.6 ± 1.3
K	15.2 ± 1.3
L	19.1 ± 1.3
M	7.6 ± 1.3
N	7.6 ± 1.3
O	20.3 ± 2.5

### Spec:

Material	CPE, 52 g
----------	-----------

### Packaging

Plastic bag	25 PCS
-------------	--------

Carton	75PCS, 3 bags, 38cm x 20cm x 31cm
--------	--------------------------------------

CE MDD 93/42/EEC

SGS

**CORE** Medical