


SMMS ISOLATION GOWN

PN: 060401 / MANUFACTURER PN: CN309

Nonwoven Isolation Gown, neck with ties, waist with ties, elastic wrists

CE MDD 93/42/EEC

SGS



Size	Universal size	
Length	A	120 ± 3cm
Width	B	140 ± 3cm
Sleeve Length	C	59 ± 2cm
Armhole Width	D	27.5 ± 2cm
Shoulder Width	E	22.5 ± 2cm
Neck Tie Length	F	35 ± 2cm
Waist Best Length	G	180 ± 5cm

Spec:

Material	SMMS, 23gsm
	EN13795 conformity

Packaging

Plastic bag	10 PCS
Carton	100PCS, 10 bags, 52cm x 30cm x 35cm