


# PP ISOLATION GOWN

PN: 060601 / MANUFACTURER PN: CN307

PP Isolation Gown, neck with ties, waist with ties, elastic wrists



CE MDD 93/42/EEC



Size	Universal size	
Length	A	115 ± 3cm
Width	B	136 ± 3cm
Sleeve Length	C	60 ± 2cm
Armhole Width	D	26 ± 2cm
Shoulder Width	E	22.5 ± 2cm
Neck Tie Length	F	35 ± 2cm
Waist Best Length	G	165 ± 5cm

### Spec:

Material	PP, 25gsm
----------	-----------

### Packaging

Plastic bag	10 PCS
Carton	100PCS, 10 bags, 52cm x 29cm x 32cm

