

PROTECTIVE GOWN

PN: 060201 / MANUFACTURER PN: CN308

PP/PE Isolation Gown, neck with velcro, waist with ties, elastic wrists

CE
0624

CE Certificate

CE 0784200148-00-00



Size	Universal size	
Length	A	120 ± 3cm
Width	B	140 ± 3cm
Armhole Width	C	27.5 ± 2cm
Sleeve Length	D	59 ± 2cm
Shoulder Width	E	22.5 ± 2cm
Waist Best Length	F	180 ± 5cm

Spec:

Material PP/PE, 42 ± 2gsm

Packaging

Plastic bag 10 PCS

Carton 100PCS, 10 bags,
52cm x 30cm x 39cm